Primo®

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SAFETY INSTRUCTIONS

- Do not clean your scale with water or immerse it into water or any other liquid.
- Never exceed the maximum weighing capacity of your scale (in such case the "OL" indication will appear on your screen). The maximum capacity of this scale is 180Kg.
- This scale is designed for household use only and not for medical or commercial use.
- This scale is a high precision measuring device. Never knock it heavily or try to disassemble it and please handle it carefully to avoid dropping or breaking during moving or during use. Step on your scale gently.
- Do not use your scale in case it malfunctions or shows any signs of damage. Take it to an authorized service center for checking, adjustment or repair.
- **Never step on your scale with wet or moist or sweaty feet. There is great danger of slipping!**
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliances by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- Use a soft slightly damp cloth in order to clean your scale and its platform. Afterwards dry thoroughly with a soft dry cloth.
- Never use hard or abrasive cleansers to clean your appliance or its platform.
- Protect your device from direct sunlight, high or extremely low temperatures, moisture and dust.
- Do not place any objects on your scale when it is off.

BATTERIES

- Insert 2 batteries type AAA by opening the battery compartment which is located on the bottom side of the scale's body and then close the battery compartment cover.
- Pay extra attention so that the batteries will be inserted with the correct polarity as shown in the battery compartment.
- Use 2 X AAA batteries of 1.5V each.
- In the case that your scale will not be used for a long time, it is advisable that you remove the batteries from your device.
- **WARNING:** Never expose the batteries in high temperatures or in direct sunlight. Do not throw the batteries in to fire. There

is danger of explosion!

- In the case that the batteries included in your scale are low and need replacement, the "LO" indication will appear on your scale.

Note:

- Set the unit of measure in "kg" or "lb" by pressing the button on the bottom side of your scale immediately after powering it on.
- The scale needs to be initialized before first use or after battery replacement. Step on the scale platform. The display will turn on, showing your weight automatically. The scale is now ready for use.

Quick weighing

If you would like to simply weigh yourself without going through the more advanced settings options, step on the scale platform, the display will turn on and show your weight automatically. And then the scale will switch off automatically after few seconds.

Using the Settings Menu for Customized Weighing

Enter your personal data (user code/height/age/gender):

1. Press the **SET** button to switch on the scale and set your personal data. In the left side of the LCD screen, "0" will begin flashing. Press the **UP (▲)** button or the **DOWN (▼)** button to select a number between 0-9 to set the data profile. Your scale can hold personal information for up to 10 different users. Press the **SET** button to confirm.
2. Next, input your height. When "165cm" starts flashing in the LCD screen, press the **UP (▲)** button or the **DOWN (▼)** button to select your height (in cm). Press the **SET** button to confirm your setting.
3. Next, input your age. When "Age 25" starts flashing in the LCD screen, press the **UP (▲)** button or the **DOWN (▼)** button to select your age. Press the **SET** button to confirm your setting.
4. Next, choose the gender of the person. Press the **UP (▲)** button or the **DOWN (▼)** button to select either male "♂" or female "♀". Press the **SET** button to confirm your gender to finish the setting.
5. When the scale shows "0.0kg", step on the scale to start measuring. Remember to remove your socks and step on the scale barefoot.
6. Stand still while the scale calculates your weight. Small square moves from the left to right on LCD screen indicating that your scale is analyzing your data. Then the scale will finish measuring, and will display weight, body fat, body water, muscle, bone,

calories, BMI, alternating between each measurement two times. Then it will automatically turn off after several seconds.

7. If no weighing is detected while "0.0" is displayed, the scale will automatically switch off after a few seconds.

Warning: Water or moist or sweat on your feet may result into a slip!

Overload indication

When the LCD displays "OL", it indicates that the maximum weight (180Kg) has been exceeded.

Attention:

- Weigh yourself without clothing and footwear, before meals and always at the same time of day.
- For better accuracy, it is not advisable to weigh yourself within 2 hours of waking.
- Keep the scale clean and dry to protect the electronic components.
- Clean after use with a slightly damp cloth, do not use solvents or immerse the product in water.
- Do not store the scale in an upright position when not in use, as this will drain the batteries.
- Avoid impact or vibration to the unit, such as dropping the scale onto the floor, or objects onto the scale.
- A user who has any metal material inserted into his/her body may get differences in measurement of body fat reading.
- **Pacemaker users should not use this device.**

Reference tables

Body fat ratio Chart (%)

AGE	Female (%)				Male (%)			
	under-weight	healthy	over-weight	obese	under-weight	healthy	over-weight	obese
6-39	<21	21~33	33~39	>39	<8	8~19	19~25	>25
40-59	<23	23~35	35~40	>40	<11	11~22	22~28	>28
60-100	<24	24~36	36~42	>42	<13	13~25	25~30	>30

Due to natural body variation between people, this chart should be used for reference purposes only.

Water Ratio Chart & Muscle Ratio Chart (%)

AGE	Female (%)				Male (%)			
	Low	Healthy	High	Body Muscle	Low	Healthy	High	Body Muscle
6-15	<21	57.0-67.0	>67.0	>34	<58.0	58.0-72.0	>72.0	>40
16-30	<23	47.0-57.0	>57.0		<53.0	53.0-67.0	>67.0	
31-60	<24	42.0-52.0	>52.0		<47.0	47.0-61.0	>61.0	
61-100	<24	37.0-47.0	>47.0		<42.0	42.0-56.0	>56.0	

Bone Weight Chart (KG)

	Female			Male		
Weight	<45 Kg	45-60 Kg	>60 Kg	<60 Kg	60-75 Kg	>75 Kg
Ref. bone weight	-2.6 Kg	-3.2 Kg	-8Kg	-3.4 Kg	-3.8 Kg	-8Kg

KCAL Chart

AGE	Female		Male	
	Weight (KG)	KCal	Weight (KG)	KCal
6-17	50	1265	60	1620
18-29	55	1298	65	1560
30-49	60	1416	70	1680
50-69	60	1302	70	1561
70-100	50	1035	65	1397

BMI Chart

BMI	Low	Healthy	High	Overweight
	<18.5	18.5-27	30-39	>40

Cleaning

The body scale must not be immersed in water or any other liquid. Wipe the scale clean using a soft, slightly dampened cloth, and then dry thoroughly using a dry soft cloth. Never use gas, hot water or other harsh chemicals to clean the scale.

Technical data

Electronic body fat scale

Model: PRBS-40464 / PRBS-40465

Battery: 2 X AAA batteries

Low battery indicator

Digital LCD display

Step on switch function

Automatic switch off

Overload indicator

Max capacity: 180 Kg / Graduation: 100g

Unit: Kg / lb

Age range for measurement: 10- 99 years

Permissible height for measurement: 100 - 240 cm

Storage spaces for different users: 10

Measuring range body fat: 5.0 – 60.0 %

Measuring range body water: 20.0 – 80.0%

Measuring range muscle mass: 25.0 – 75.0%

Measuring range bone mass: 1.0 – 10.0%

Measuring range Calories: 500 - 5000 Kcal

Weighing scope: 5Kg to 180Kg max

This product is in strict conformity with all the valid directives of European Union applying on this type of electrical appliances.

Instructions for the correct disposal of the unit according to European directive 2002/96/EU

Should you find one day that your appliance needs to be replaced or if it is of no further use to you, think of the protection of the environment:

- 1) Do not dispose your appliance along with the rest of the public waste (this is also the meaning of the shown recycling sign).
- 2) Contact your Public Authorities and they will instruct you of the recycling centers to which your appliances must be disposed.
- 3) Correct disposal of your appliance helps the protection of the environment as well as the recycling of the appliance components.
- 4) The electrical appliances due to their construction materials, if not correctly disposed, may lead to environmental and furthermore health hazards.

BATTERIES DISPOSAL

Used batteries should not be disposed along with the rest of the domestic waste. Dispose of your batteries into the battery collection bins or contact your Public Authorities for advice on the battery collection points.

With the warranty of:

Buy Way

BUY WAY S.A.
Tyrtaiou Avenue,
31st km Athinon-Lamias,
19014 Afidnes, Greece
Tel.: +30 210-2464214
www.buyway.gr